

# Vented 'Cirri' Reusable <u>Full & Nasal Face</u> CPAP Mask with Headstrap



### Vented Reusable 'Cirri' Full Face CPAP Masks with Headstrap (PC)(S)

- 3D mask cushion provides an excellent seal without discomfort
- Cushion is easily attached and detached from the frame
- Quick release clips for easy removal of the headgear
- Innovative adjustable forehead support provides one handed adjustment and stability
- Flexible, one piece forehead pads evenly conform to the contour of the forehead
- Quiet air vent minimises noise and enables efficient CO<sub>2</sub> washout and reduces patient disturbance while venting air

### **Part Numbers**

PH11343 Vented Reusable 'Cirri' Full Face CPAP Masks with Headstrap - **Small**PH11344 Vented Reusable 'Cirri' Full Face CPAP Masks with Headstrap - **Medium**PH11345 Vented Reusable 'Cirri' Full Face CPAP Masks with Headstrap - **Large** 



### Vented Reusable 'Cirri' Nasal Face CPAP Masks with Headstrap (PC)(S)

- 3D mask cushion provides an excellent seal without discomfort
- Cushion is easily attached and detached from the frame
- Quick release clips for easy removal of the headgear
- Innovative adjustable forehead support provides one handed adjustment and stability
- Flexible, one piece forehead pads evenly conform to the contour of the forehead
- Quiet air vent minimises noise and enables efficient CO<sub>2</sub> washout and reduces patient disturbance while venting air

### **Part Numbers**

PH11353 Vented Reusable 'Cirri' Nasal Face CPAP Masks with Headstrap - **Small**PH11354 Vented Reusable 'Cirri' Nasal Face CPAP Masks with Headstrap - **Medium**PH11355 Vented Reusable 'Cirri' Nasal Face CPAP Masks with Headstrap - **Large** 





# Non-Vented Reusable <u>Full & Nasal Face</u> CPAP Mask with Headstrap



# Non-Vented Reusable Full Face CPAP Masks (S) (PC)

- Suitable for patients suffering from acute or chronic respiratory distress or failure
- Specifically designed for use with ventilators that use an active exhaust valve

#### **Part Numbers**

PH10614 Non-vented Reusable 'Ease' Full Face Masks - **Small**PH10615 Non-vented Reusable 'Ease' Full Face Masks - **Medium**PH10616 Non-vented Reusable 'Ease' Full Face Masks - **Large** 



## Non-Vented Reusable Nasal Face CPAP Masks (S) (PC)

- Suitable for patients suffering from acute or chronic respiratory distress or failure
- Specifically designed for use with ventilators that use an active exhaust valve

### **Part Numbers**

PH10604 Non-vented Reusable 'Ease' Nasal Face Masks - **Small**PH10605 Non-vented Reusable 'Ease' Nasal Face Masks - **Medium**PH10606 Non-vented Reusable 'Ease' Nasal Face Masks - **Large** 





