






AUSTRALIAN OWNED & OPERATED







Company Profile

Founded in 1985, Parker Healthcare Pty Ltd is an Australian owned and operated medical supply company, with offices located in all main states of Australia and international offices in both New Zealand and Malaysia.

As a Quality Assured ISO 9001:2015 accredited company, Parker Healthcare supplies a diverse range of products that covers the neonatal through to adult divisions of the hospital. These segments include, but are not limited to; Anaesthesia and Respiratory Care, Infant and Paediatric Care, Maternity, Emergency and ICU.

Parker Healthcare have specialised divisions covering Hospital, Pharmacy, Aged Care, Veterinary and Home Healthcare, enabling the purchase and supply of products specific to these industries.

Parker Healthcare is committed to working towards continual improvement in the quality of care given to all patients. Through understanding our customers' needs and expectations this enables us to develop consumables and equipment to address the changing needs of both Australian and International healthcare.

Core Capabilities

- Own manufacturing facility in Malaysia (approx. 32,000 sq ft)
- Outstanding customer service
- Wide range of products
- · Stringent quality control checks
- Highly experienced and trained workforce with over 35 years of industry experience
- Service Department for fast and local service & maintenance of goods
- ISO 9001:2015 Accredited

Mr. Ian Bird Managing Director















Contents










Page 16



























VISIOFOCUS® Non-contact Thermometer with projected temperature image

The only non-contact thermometer in the World that takes temperature and projects it onto the surface being scanned (including liquids)

- Environment and ambient working temperature: 16 40°C
- Forehead measuring range: 34 42.5°C
- General measuring range: 1 55°C
- Quick calibration with MQCS & AQCS
- Projected image indication of correct reading distance.
- Silent, non-contact reader does not disturb your child's sleep.
- Measures and projects the temperature of objects, food, liquids and the skin, as well as the ambient temperature.
- 5 Colour illuminated display.
- 9 Recallable memories.

Part No Description

06400 Visiofocus Non-contact Image Projector Thermometer

TGA #125945



















BodyTemperature

Rottle

Food

Rath

Bedding

Skin Temperature

Ambient

Memory Recall



VISIOFOCUS® Mini Non-contact Thermometer with projected temperature image

The only non-contact thermometer in the World that takes temperature and projects it onto the surface being scanned (including liquids)

- · Projected image indication of correct reading distance.
- · Silent, non-contact reader does not disturb your child's sleep.
- Measures and projects the temperature of objects, food, liquids and the skin, as well as the ambient temperature.
- Forehead measuring range: 34 42.5°C
- General measuring range: 1 80°C

Part No Description

06700 Non-contact Image Projector Thermometer

TGA #125945





















Thermometers



PARKERmed Non-contact Infrared Thermometer

- · Comparatively lightweight and small in size
- · Adopting reliable sensor made in Germany
- · All-new & Patented probe design to ensure high accuracy
- Measurement mode for human body and object available. One button to switch measurement mode
- Beeper for fever or high temperature
- · 3 Colours Backlight LCD for easy reading
- · C / F reading available
- 10 memories
- · Automatic shut off

Part No Description

PM9000 Non-contact Infrared Thermometer

TGA #336632














PM1000 Non-Contact Infrared Thermometer

- · Comparatively lightweight and small in size
- Rapid Temperature Measurement
- · Fever Alert Setting
- Stores Last 10 Readings
- Motor Vibration Reminder
- · High Quality Sensor for Accurate Temperature Measurement
- Two Measurement Modes for Body or Object
- · Automatic Shut off
- · Quality LCD Screen


Part No Description

PM1000 Non-contact Infrared Thermometer

TGA #336632









Silicone Fixation Tape





Can also be used to hold down face mask to stop glasses fogging up See the results below







RESULTS WITHParker Sil-Flex™ Silicone Tape

Parker Sil-Flex™ Silicone Fixation Tape

Applications:

- · Fragile / Sensitive Skin: Neonates to Geriatric patients.
- Suitable on damaged skin such as Dialysis, Chemotherapy and Radiography patients.
- · All types of fixation e.g. IV, Cannula, Tubes, Electrodes etc.

Product Features:

- · Excellent adhesive property, leaves no residue.
- · Does not stick to hair.
- · Waterproof.
- · Suitable for Fragile / Sensitive Skin.
- · Can be re-applied without losing the adhesiveness.
- Patented design, allowing you to tear to the desired length (perforated every 1cm) without compromising the tensile strength of the tape (no scissors required).
 (Applies only to PH-ST1519, PH-ST1525 & PH-ST5025)
- Transparent, soft and pliable, will not damage or disrupt any skin layers.
- Higher conforming property due to PU backing when compared to existing tapes.

| Part No. PH-ST1519 | Description 1.5m Length Roll, 1.9cm Width, Perforated | QTY Box of 6 rolls |
|--------------------|---|-----------------------|
| PH-ST1525 | 1.5m Length Roll, 2.5cm Width, Perforated | Box of 6 rolls |
| PH-ST5025 | 5.0m Length Roll, 2.5cm Width, Perforated | Box of 6 rolls |

PH-STS2050 5cm Length, 2.0cm Width, 2 strips/pack Box of 50 packs Beige



PH-ST1519 1.5m Length, 1.9cm Width



PH-ST1525 1.5m Length, 2.5cm Width



PH-ST5025 5.0m Length, 2.5cm Width



PH-STS2050 5cm Length, 2cm Width. Beige





Blood Pressure Monitor & Blood Glucose Meter









Reflects Patient Blood Pressure in color


PARKERmed Talking Blood Pressure Monitor

- · Large LCD Display with Tricolour backlight
- · Memory function (2 users x 120 groups of measured values)
- Alarm function
- · Auto power off 3 mins after last operation
- Irregular heart beat detection function
- Easy to operate, automatic operation
- LCD Display emits colour to reflect Low, Normal or High Blood Pressure reading

Part No Description

PM-BP386A Talking Blood Pressure Monitor





Blood Glucose Meter

- CE0123, IS013485, GMP
- · No Coding
- · OEM Welcomed
- Accurate, Quick, Easy-to-use Blood Glucometer
- · Alternative site testing
- · Wide measuring range of 20 to 600 mg/dL
- Fast testing time in 6 seconds
- Only 0.7µL blood volume
- 120 test memory with date and time
- 7, 14 and 28 days average

Part No Description

GM-01 Blood Glucose Meter - Complete Set

GM-01A Blood Glucose Test Strip (Package of 25pcs strips)





| Range of Glucose Measurement | 20-600 mg/dL (1.1-33.3mmol/L) |
|-----------------------------------|--|
| Sample Type | Capillary Whole Blood |
| Calibration against (sample type) | Plasma - Equivalent |
| Test Time | 6 seconds |
| Sample Volume | 0.7μL |
| Sample Application | Capillary Effect |
| Operating Temperature | 10-50°C (50~104°F) |
| Operating Humidity | 20-80% RH |
| Memory Capacity | 180 |
| Risk Indicator | SMILE when test within 70-120mg/dL (3.9~6.6mmol/L) POUT when test <70 mg/dL (3.9mmol/L) or >180mg/dL (10mmol/L) |
| Average | 7,14 and 28 days |
| Downloadable | RS-232 |
| Battery Type | V Li-Battery (CR2032) x 1 |
| Battery Life | More than 1,000 tests |
| Auto OFF | 3 minutes |
| Meter Dimensions | 80 x 48 x 15mm |
| Meter Weight | 37g (including battery) |



Finger Pulse Oximeters





Bluetooth Adult Finger Pulse Oximeter

- · Highly Accurate
- Display Sp02, PR, PI
- · Colour OLED Display with 2 modes
- · Suitable for Adult, Adolescent, Child & Infant
- · Automatic power-off design
- · Protective Carry Case

Part No Description

PH-MD300C1218

Adult Finger Pulse Oximeter

TGA #291638



Oxywatch Adult SpO_2 Finger Pulse Oximeter

- Dual colour OLED displays SpO2, PR, Pulse bar, and waveform
- · Low power consumption; real-time battery status indication
- Level 1-10 adjustable brightness.
- 6 display modes
- · 2pcs AAA-size batteries; battery-low indicator
- · Automatically power off

Part No Description

PH-MD300C2F Adult Finger Pulse Oximeter



Adult SpO₂ Finger Pulse Oximeter

- Highly Accurate
- Display Sp02 waveform (pleth or line), Sp02, PR, bar graph
- · LCD display screen with six display modes
- Automatic power-off design

Part No Description

PH-MD300C22 Adult Finger Pulse Oximeter

TGA #291638



Adult/Paediatric SpO₂ Finger Pulse Oximeter

- · Highly Accurate
- Display Sp02 waveform (pleth or line), Sp02, PR, bar graph
- · LCD display screen with six display modes
- · Automatic power-off design

Part No Description

PH-MD300-C5 Adult/Paediatric Finger Pulse Oximeter

TGA #291638



Paediatric SpO₂ Finger Pulse Oximeter

- Highly Accurate
- · Special cartoon design for children
- Display Sp02 waveform (pleth or line), Sp02, PR, bar graph
- LCD display screen with six display modes
- · Automatic power-off design

Part No Description

PH-MD300C53 Paediatric Finger Pulse Oximeter (Bear)
PH-MD300C53 Paediatric Finger Pulse Oximeter (Frog)
PH-MD300C55 Paediatric Finger Pulse Oximeter (Panda)

TGA #291638



Nebulizer





Portable Nebulizer

- Vibrating Mesh provides fine mist
- 5 µm particale spray above 60%
- · Compact and Portable
- Automatic shutdown after 20 minutes
- Built-in Low Battery indicator
- · Power options: 2 AA batteries or USB cable

Part No Description

PH-NBM-1 Portable Mesh Nebulizer

Spare Parts:

Part No Description **NML-006** Large mask NMS-004 Small mask NMP-002 Mouthpiece

PH-NBM-001 **USB** Power adaptor









Tiny Particle

Nebulizer

Reliable Quality

Easy to Use

Specifications:

| Power Consumption | < 1.5W |
|-------------------------------------|--|
| Working current consumption | < 500 mA |
| Sound level | Less than 50 dB at 1 metre |
| Operation environment | 10 ~ 40°C (50°F ~ 104°F), 30% ~ 80% RH |
| Storage environment | -20 ~ 55°C (-4°F ~ 131°F). ≤85% RH |
| Nebulizers rate | ≤ 0.2mL/min |
| Time AUTO-OFF | Power on, operation 20 min to auto-off |
| Medication bottle volume | 6 mL |
| Ultrasonic Oscillation Frequency | 103 KHZ ~ 120 KHZ ± 3K |
| Particle Size | 5μm or less accounted for more than 60% |
| Battery life | Continuous operation ≥ 1 hour, (1.5V X2, AA, Alkaline) |
| Standards | EN60601-1-1, EN60601-1-2 |

Accessories:



Large mask

Part No: NML-006



Small mask

Part No: NMS-004



Mouthpiece



USB Power Adaptor

Part No: NMP-002





Oxygen Concentrator PARKER







Zen-O™ can hold up to two 12 cell batteries



Zen-O™ with bag and pull cart

Zen-O™ Portable Oxygen Concentrator

- Dual FlowMode enables patient to alternative between continuous flow and pulse mode oxygen therapy
- Simple to use Intuitive button operation
- LCD Display
- · Rate Responsive Therapy oxygen delivered on demand
- triggered by the patient's inhalation
- Eco Mode switch between Rate Responsive to Fixed Minute
- · Volume for longer battery duration
- Auto Mode activates after 60 seconds if no breath is detected and automatically deliver pulses at a rate of 18 breaths per minute
- · Easily replaceable sieve modules
- Durable and reliable

RS-00523

| Part No | Description |
|----------------|---|
| RS-00502-AUS-D | Zen-O™ concentrator - 2 battery package |
| RS-00501 | Zen-O™ battery |
| RS-00509 | Zen-O™ carry bag |
| RS-00507 | Zen-O™ cart |
| RS-00508 | Zen-O™ DC adaptor |
| RS-00511 | POC filter wrench |
| RS-00512 | POC cannula filter (packet of 10) |
| RS-00513 | Sieve module assembly |
| RS-00520 | Zen-O™ AC power supply |

Each Zen-O $^{\text{TM}}$ device (Part No: RS-00502-AUS-D) is supplied with a carry bag, two batteries, AC/DC power cables and a pull cart.

Accessories Bag



Accessories bag



Power cords, AC and DC power supply









Zen-O™ lite carry bag



Power cords, AC and DC power supply

Zen-0 /ite™ Portable Oxygen Concentrator

- · Compact, portable and lightweight weighing just 2.5kg
- Simple to use Intuitive button operation
- LCD Display
- Rate Responsive Therapy oxygen delivered on demand triggered by the patient's inhalation
- Eco Mode switch between Rate Responsive to Fixed Minute Volume for longer battery duration
- Auto Mode activates after 60 seconds if no breath is detected and automatically deliver pulses at a rate of 18 breaths per minute
- Easily replaceable sieve modules
- · Durable and reliable

| Part No | Description |
|----------------|--|
| RS-00608-AUS-D | Zen-O lite™ concentrator with two 8 cell batteries |
| RS-00601 | Zen-O lite™ rechargable battery |
| RS-00602 | Zen-O lite™ AC power supply w/EU cord |
| RS-00605 | Zen-O lite™ Dc power supply |
| RS-00606 | Zen-O lite™ carry bag |
| RS-00616 | Replacement sieve modules |
| RS-00617 | Zen-O lite™ cannula wrench |
| RS-00619 | Zen-O lite ™ backpack |
| RS-00512 | Zen-O lite™ cannula filter (packet of 10) |
| RS-00523 | Accessories Bag |

Each Zen-O lite $^{\text{TM}}$ device (Part No: RS-00608-AUS-D) is supplied with a carry bag, two batteries, AC/DC power cables.



CPAP Masks & Smooth-Bor Tubing






"Vented" Reusable 'Cirri' Full Face CPAP Masks with Headstrap

- NEW 3D mask cushion provides an excellent seal without tight headgear and discomfort
- · Simple action to attach and detach the cushion from frame
- · Quick release clips allow for easy removal of the headgear
- Innovative adjustable forehead support provides one-handed adjustment and stability
- Flexible, one-piece forehead pads evenly conform to the contour of the forehead
- Quiet air vent that minimizes noise and enable efficient CO2 washout and vented air is spread as a cloud so that it does not disturb patient or companion


| Part No PH11343 | <i>Description</i> Full Face | <i>Size</i> Small | <i>QTY</i> Each |
|------------------------|------------------------------|----------------------|--------------------|
| PH11344 | Full Face | Medium | Each |
| PH11345 | Full Face | Large | Each |



"Vented" Reusable 'Cirri' Nasal Face CPAP Masks with Headstrap

- NEW 3D mask cushion provides an excellent seal without the discomfort of tight headgear
- · Simple action to attach and detach the cushion from frame
- · Quick release clips allow for easy removal of the headgear
- Innovative adjustable forehead support provides Innovative stability
- Flexible, one-piece forehead pads evenly conform to the contour of the forehead
- Quiet air vent that minimizes noise and enable efficient CO2 washout and vented air is spread as a cloud so that it does not disturb patient or companion.

| Part No | Description | Size | QTY |
|---------|-------------|--------|------|
| PH11353 | Nasal Face | Small | Each |
| PH11354 | Nasal Face | Medium | Each |
| PH11355 | Nasal Face | Large | Each |



Disposable Adult Smooth-Bor Tubing - 22mm

| Part No | Colour | Length |
|-------------|---------|-------------|
| MP-SB22N48 | Natural | 1.2m (48") |
| MP-SB22N72 | Natural | 1.8m (72") |
| MP-SB22N96 | Natural | 2.4m (96") |
| MP-SB22N120 | Natural | 3.0m (120") |
| | | |

*Other styles and lengths available upon request. Suitable for use with CPAP Machines.



Reusable Adult Smooth-Bor Tubing with Cuff - 22mm

| Part No | Length |
|-----------|------------|
| SB403-48 | 1.2m (48") |
| SB403-72 | 1.8m (72") |
| SB403-96 | 2.4m (96") |
| SB403-120 | 3.0m (120" |


^{*}More lengths available upon request. Suitable for use with CPAP Machines.











Preemie Muffs™ **Premature Neonatal Ear Covers**

Protecting our smallest patients from acoustic noise in the NICU.

It is now known that loud noises can negatively impact on physiological parameters such as blood pressure, breathing, heart beat and oxygen saturation. When exposed to loud noises, a baby's sleep is interrupted, which then increases stress and energy consumption, which means less energy is used for growth and development.

In the neonatal intensive care unit and during transportation, a premature baby is exposed to a variety of noises. The Parker Neocare Preemie Muffs™ are a safe and comfortable way to provide hearing protection from noises that are generated from:

- · Ventilation machines or open care systems
- · External and Internal transfers
- · During procedures, such as cardiac echoes
- · Acoustic noises from the Neonatal Nursery when incubator doors are open or when the neonate or infant is outside the incubator.

Preemie Muffs™ are the safe and comfortable way to provide hearing protection for premature babies. They are commonly used within loud environments in the neonatal nursery and during transportation.

Part No. Description QTY **MP-EM200** Preemie Muffs™ Ear Covers 25 pairs/box



Nasal Cannulas





Oxygen Nasal Curved Tips Cannula

- · Contoured soft cannula provides total patient comfort
- · Star lumen tubing
- · Over the ear fit secures position on the patient
- Supplied as box of 50

| Part No. | Size | lube Length | QIY |
|--------------|------------|------------------|-----|
| PH2801 | Adult | 2.1 metre | 5 |
| PH2801-GR | Adult | 2.1 metre, Green | 5 |
| PH2801-3GR | Adult | 3 metre, Green | 5 |
| PH2803 | Paediatric | 2.1 metre | 5 |
| PH2803-3GR | Paediatric | 3 metre, Green | 5 |
| PH2804 | Infant | 2.1 metre | 5 |
| PH2804-3GR | Infant | 3 metre, Green | 5 |
| PH2805 | Neonate | 2.1 metre | 5 |
| PH2805-3GR | Neonate | 3 metre, Green | 5 |
| Flared Tip | | | |
| PH2801-F | Adult | 2.1 metre | 5 |
| PH2801-F-GR | Adult | 2.1 metre, Green | 5 |
| PH2801-F-3GR | Adult | 3 metre, Green | 5 |



"Softi" Smoothflow® Oxygen Nasal Curved Tips Cannula

- Contoured "Softi" Smoothfl ow® cannula provides total patient comfort
- · Over the ear fit secures position on the patient

| Part No | Size | Connection | QTY |
|---------|------------|---------------------|-----|
| DR-0161 | Adult | 6mm M Barbed Nipple | 5 |
| DR-0163 | Paediatric | 6mm M Barbed Nipple | 5 |



"Softi" Smoothflow® Diagnostic Cannula, Softi Smoothflow, Oxygen, Nasal Curved Tips Cannula

- Can be used in conjunction with pressure transducers or thermistors in Polysomnography (Sleep Studies)
- Optimizes pressure/snore waveform signals
- Inexpensive
- Contoured "Softi" Smoothflow® cannula provides total patient comfort
- Over the ear fit secures position on the patient
- · Available in various lengths
- Male luer lock connection

| Part No. | Size | Length | Connection | QTY |
|-----------|------------|--------|------------|-----|
| LD50ML61 | Adult | 50cm | Male Luer | 5 |
| LD50ML63 | Paediatric | 50cm | Male Luer | 5 |
| LD100ML61 | Adult | 100cm | Male Luer | 5 |
| LD100ML63 | Paediatric | 100cm | Male Luer | 5 |
| LD130ML61 | Adult | 130cm | Male Luer | 5 |
| LD130ML63 | Paediatric | 130cm | Male Luer | 5 |

^{*}Available with Female Luer upon request



Oxygen Tubing



Non-Kink Star Lumen Disposable Oxygen Tubing

- · Kink resistant tubing with inner star lumen design
- Flexible PVC tubing with standard connectors
- Lengths available in 2.1, 3, 4, 6, 10, 15 and 30 metre variations

| Description | QTY |
|---|--|
| 2.1 metre oxygen tubing | 5 |
| 2.1 metre oxygen tubing GREEN | 5 |
| 3 metre oxygen tubing | 5 |
| 3 metre oxygen tubing GREEN | 5 |
| 4 metre oxygen tubing | 5 |
| 6 metre oxygen tubing | 5 |
| 10 metre oxygen tubing GREEN | 2 |
| 15 metre oxygen tubing GREEN | 2 |
| 10.0m, 6mmM, 4mmF Green Oxygen Tube, White Cuffs | 2 |
| 15.0m, 6mmM, 4mmF Green Oxygen Tube, White Cuffs | 2 |
| | 2.1 metre oxygen tubing 2.1 metre oxygen tubing GREEN 3 metre oxygen tubing 3 metre oxygen tubing GREEN 4 metre oxygen tubing 6 metre oxygen tubing 10 metre oxygen tubing GREEN 15 metre oxygen tubing GREEN 10.0m, 6mmM, 4mmF Green Oxygen Tube, White Cuffs 15.0m, 6mmM, 4mmF |



Flexible Oxygen Bubble Tubing

- · Flexible oxygen bubble tubing
- · Available in 30 metre rolls
- Can be cut to any length as required
- · Fits multiple oxygen therapy devices and connectors

OBTG-3-30 Oxygen Bubble Tubing

3.0 mm ID, Bubble every 90cm,

(30m per roll)



Face Protective Shield - Anti Fog

- Protects facial area and associated mucous membranes (eyes, nose, mouth) from splashes and sprays of body fluids.
- Light Weight, Durable and Comfortable
- Suitable for use by healthcare workers, dental providers, veterinary care
 personnel, laboratory workers, pre-hospital emergency medical providers,
 police, firefighters, and custodial staff dealing with germs, spills and
 contaminated waste.

| Part No | Description | QTY |
|------------|---------------------------------|-----|
| FPS02-AF-P | Anti Fog Face Protective Shield | 5 |







Reusable Neonatal / Infant Silicone Donut Masks

Part NoSizeUOMPH10111Size 1, Large Infant, 15mm 0DEachPH10112Size 2, Small Paediatric, 22mm IDEach

Reusable and Autoclavable up to 134°C



Oxygen Therapy Paediatric <u>Elongated</u> Mask with Oxygen Tubing

Disposable <u>Elongated</u> O2 Therapy Mask with Elastic Strap and 2.1m Oxygen Tubing Mask Colour: Clear/Transparent

Part No Description QTY

PH0166-2.1 <u>Elongated Paediatric (2.1m Clear Tubing)</u> Each

PH0166-2.1GR Elongated Paediatric (2.1m Green Tubing) Each



Oxygen Therapy Adult Elongated Mask with Oxygen Tubing

Disposable <u>Elongated</u> O2 Therapy Mask with Elastic Strap and 2.1m Oxygen Tubing Mask Colour: Clear/Transparent

Part No Description QTY

PH0188-2.1 Elongated Adult (2.1m Clear Tubing) Each

PH0188-2.1GR Elongated Adult (2.1m Green Tubing) Each



Aerosol Therapy Paediatric & Adult Masks With Tubing

Disposable Aerosol Therapy Mask with Elastic Strap, 10ml Nebulizer Bottle, 2.1m Oxygen Tubing Mask Colour: Clear/Transparent

| Part No | Description | QTY |
|-------------|--|------|
| PH0266N-2.1 | Elongated Paediatric (2.1m Clear Tubing) | Each |
| PH0266N-3GR | Elongated Paediatric (3m Green Tubing) | Each |
| PH0288N-2.1 | Elongated Adult (2.1m Clear Tubing) | Each |
| PH0288N-3GR | Elongated Adult (3m Green Tubing) | Each |
| | | |



Tracheostomy Masks, Filters & Nebuliser Bottle



Tracheostomy Masks with 6mm Connector

Part No Description QTY

Each

PH30162 Disposable Adult O₂ Tracheostomy Mask with Elastic Strap & 6mm Connector

PH30163 Disposable Paediatric O₂ Tracheostomy Mask Each

with Elastic Strap & 6mm Connector

HME Trach





| Description | HME Trach T | HME Trach Basic |
|-----------------|--------------------|--------------------|
| Part Number | PH51050902 | PH51050901 |
| Connector | 15F | 15F |
| Dead Space | 8ml | 14ml |
| Pressure Drop | 0.15KPa@30L/Min | 0.15KPa@30L/min |
| Moisture output | ≥26mg/L @ VT=500ml | ≥26mg/L @ VT=500ml |
| Weight | 7g | 9g |
| Box quantity | 100 | 100 |



Nebuliser Bottle

Part No QTYDescription PH-NBG18 Nebuliser Bottle, 18mm OD, Each


6ml Volume



Resuscitators & Connectors



Each



Part Number: PH30038

Description: Swivel Oxygen Connector

Accepts 3mm to 7mm ID Tubing

UOM: Box Qty:



Part Number: MP58003

Description: Multi-size Tubing Joiner

Accepts 5mm to 8mm ID Tubing

UOM: Each Box Qty: 5





With POP-OFF Valve

Disposable Manual Resuscitator with POP-OFF Valve

Part No Description QTY

Adult, 1500ml Ventilation Bag, 60cmH $_2$ O Pop-off Valve, Size 5 Disposable Mask, 2500ml Reservoir Bag and 2.1m Non-Kink Tubing PH60101

PH60102 Paediatric, 550ml Ventilation Bag, 40cmH₂O Pop-off Valve, Fach

Size 3 Disposable Mask, 2500ml Reservoir Bag

and 2.1m Non-Kink Tubing

PH60103 Infant, 280ml Ventilation Bag, 40cmH₂O Pop-off Valve, Each

Size 1 Disposable Mask, 600ml Reservoir Bag

and 2.1m Non-Kink Tubing



CPR Resuscitators

Part No Description

PH30100 CPR Resuscitator with Case

(Mask, one way valve and filter in a plastic case)

PH30103 Pocket CPR Resuscitator with Case on a Keyring

(Shield and Valve/Filter in a plastic pocket case)





Roll-on Pain Relief Cream

Convenient, easy to apply roll-on pain relief for sprains, swelling, backache and headache. Also relieves itchy skin and stimulates blood circulation.

- · Herbal extract
- · Quick pain relief
- · Sanitary & convenient ball design
- · Hands free application
- · Easily absorbed gel
- · Airtight container
- Contains Menthol, 4% Lidocaine, Wintergreen Oil, Eucalyptus Oil, Camphor, Clove Oil, Turpentine Oil, Thermal agent.

Part No Description

PH-QS8104-H Pain Relieving Gel Roll-On

Heavy (Extra Strength)

PARKER

PAIN RELIEVING

ROLL-ON

Extra Strength

ROLL-ON

Colour: Blue

PH-QS8104-L Pain Relieving Gel Roll-On

Light Strength Colour: Red





Insect Repellent & Relief





Rotating Mosquito Repellent Bracelet

Attractive child-friendly rotating bracelet containing natural plant based mosquito repellent, infused with essential oils protecting infants & children up to 1.5 metres.

- Easy to wear Silica Gel Snap Bracelet
- Waterproof & Sweat Proof PVC, PP, ABS material
- · Removable bracelet, attachable to collar or backpack

Part No Description



PH-QS9007-Y Rotating Mosquito Repellent

Bracelet - Yellow



PH-QS9007-P Rotating Mosquito Repellent

Bracelet - Pink



PH-QS9007-B Rotating Mosquito Repellent

Bracelet - Blue



















Mini-Click® Insect Bite Relief

The patented safe, quick and easy way to stop that urge to scratch. Mini-Click® is a revolutionary and effective method for eliminating irritation of insect bites. Using electric micro-impulses to reduce histamines locally, irritation is rendered inactive by the application of electric pressure.

- · Safe for adults and children over 4 years
- · Clinically tested
- Non allergic
- · Italian designed and patented
- Functions with piezoelectricity, no batteries required

Part No Description
PH-06300 Mini-Click®

Maximum output voltage: 2.24kV

Maximum output current under a load resistance of 10,000ohm: 199.7mA

Duration of single discharge: 0,01µs
Dimensions: 24mm x 18.8mm x 42.4mm

Weight: 5 grams

Piezoelectricity: 20,000 discharges

Blister Pack Dimensions: 80mm x 180mm Weight: 11g



Neonatal First Aid Kit & Portable Toilet



Neonatal First Aid Kit

A first aid kit for basic wound treatment in emergency situations.

Application

For emergency treatment of accidental injuries

Instructions

Open the kit, clean the wound, choose the suitable bandage or dressing to cover the wound. Fix the dressing with a bandage or tape if necessary.

Contraindications

Patients who are allergic to any of the components should not use, if allergic reaction occurs, stop use immediately and consult a doctor if necessary.

Products are all single use.

Part No. Description

FK-023 Neonatal First Aid Kit

Products Included:

| Item | Size | Qty |
|------------------------------|---------------|--------|
| PE Tape | 1.25cm x 9.1m | 1 roll |
| PET Bandage | 5 x 250cm | 1 roll |
| Alcohol Pad | 3 x 6cm | 6 pcs |
| Hydrocolloid Blister Plaster | 41 x 62mm | 1 pc |
| Hydrocolloid Blister Plaster | 20 x 60mm | 2 pcs |
| Band-aids (Elastic Fabric) | 76 x 38mm | 1 pc |
| Band-aids (Elastic Fabric) | 72 x 19mm | 4 pcs |
| Band-aids Waterproof | 72 x 19mm | 4 pcs |



Baby and Kid's Portable Toilet



The container can be extended (400ml)

Happy-Pi Portable Toilet

Flexible and extendable body for easy storage

- Unisex
- Hygenic
- · Comes with a colourful design
- Three character faces: Cat, Froggy and Ducky

PM-MI171-Y Portable Toilet - Yellow

PM-MI171-W Portable Toilet - White

Ear Wax Remover & Security Ties





Safe & flexible suction nozzle






Ear Wax Remover Vacuum Cleaner Kit

The Ear Wax Remover Vacuum Cleaner is perfect for removing ear wax, gently vibrating and vacuuming wax, dirt particles and moisture from the ears.

- · Safe and simple to use
- · Painless and super effective
- Removes wax without injury
- Vibrates and vacuums simultaneously
- Durable, lightweight and long-lasting
- Convenient and portable
- Suitable for adults & children
- High quality material

Part No Description

PH-ED100 Ear Wax Remover Vacuum Cleaner Kit





















SECUR-A-TAG Security Ties

- · Tamper-Evident Design
- · Assortment of Colours
- · Easy to Use Loop, Insert & Pull
- · Write-on Identification Tag
- · Secures up to 45mm diameter
- · Hundreds of Uses

73287-50/A Security Ties (Assorted Colours)

10 of each colour (Box of 50)

73287-50/B Security Ties

BLUE only (Box of 50)

73287-50/G Security Ties

GREEN only (Box of 50)

73287-50/P Security Ties

PINK only (Box of 50)

73287-50/R Security Ties

RED only (Box of 50)

73287-50/Y Security Ties

YELLOW only (Box of 50)

NOTE: Colours may vary slightly from pictured items





SECUR-A-TAG The essential travelling companion



Index

| В | | Neonatal First Aid Kit | 21 |
|--|----|---|-------|
| Blood Glucose Meter | 8 | Non-Kink Star Lumen Disposable Oxygen Tubing | 16 |
| Blood Pressure Monitor | 8 | 0 | |
| C | | Oxygen Concentrator, Portable | 11–12 |
| Cannula, Diagnostic, Softi Smoothflow, Oxygen, Nasal Curved Tips | 15 | Oxygen Therapy Masks | 17 |
| Cannula, Oxygen, Nasal Curved Tips | 15 | Oxygen Tubing | 16 |
| Cannula, Softi Smoothflow, Oxygen, Nasal Curved Tips | 15 | P | |
| Connectors/Joiners | 19 | Paediatric SpO2 Finger Pulse Oximeter | 9 |
| CPAP Masks | 13 | Pain Relief Cream, Roll On | 19 |
| CPR Resuscitators | 18 | Parker Sil-Flex [™] Silicone Fixation Tape | 7 |
| D | | PARKERmed Non-contact Infrared Thermometer | 6 |
| Disposable Oxygen Tubing | 16 | PARKERmed Talking Blood Pressure Monitor | 8 |
| E | | Portable Nebulizer | 10 |
| | 22 | Portable Oxygen Concentrator | 11–12 |
| Ear Wax Remover, Vacuum Cleaner Kit | 22 | Portable Toilet, Happy-Pi | 21 |
| F | | Preemie Muffs, Ear Covers, Neonatal | 14 |
| Face Protective Shield | 16 | Protective Face Shield | 16 |
| Finger Pulse Oximeters | 9 | Pulse Oximeters, Finger, Adult/Paediatric | 9 |
| First Aid Kit | 21 | R | |
| Fixation Tape, Silicone | 7 | Resuscitator, Manual, Disposable | 19 |
| Flexible Oxygen Bubble Tubing | 16 | Resuscitators, CPR | 19 |
| Н | | Reusable Smooth-Bor Tubing | 13 |
| HME Trach Filters | 18 | S | |
| I control of the cont | | Security Ties, Secur-A-Tag | 22 |
| Insect Bite Relief, Mini-Click | 20 | Silicone Fixation Tape | 7 |
| M | | т | |
| Masks, Aerosol Therapy, Disposable | 17 | Talking Blood Pressure Monitor | 8 |
| Masks, Cirri, CPAP, Full Face with Headstrap, Reusable | 13 | Tape, Silicone, Fixation | 7 |
| Masks, Donut, Silicone, Neonatal/Infant, Reusable | 17 | Thermometers | 5–6 |
| Masks, Oxygen Therapy, Elongated, Disposable | 17 | Toilet, Happy-Pi Portable Toilet | 21 |
| Masks, Oxygen Therapy, Elongated, Reusable | 17 | Tubing, Oxygen, Bubble, Flexible, Disposable | 16 |
| Masks, Tracheostomy, Reusable | 18 | Tubing, Oxygen, Non-Kink Star Lumen, Disposable | 16 |
| Masks, Vented, Cirri, CPAP, Full Face with Headstrap, Reusable | 13 | Tubing, Smooth-Bor, Adult, Disposable | 13 |
| Mosquito Repellent Bracelet | 20 | Tubing, Smooth-Bor, Adult, Reusable | 13 |
| N | | V | |
| Nasal Face CPAP Masks | 13 | VISIOFOCUS® Mini Non-contact Thermometer | 5 |
| Nebuliser Bottle | 18 | VISIOFOCUS® Non-contact Thermometer | 5 |
| Nebulizer, Portable | 10 | Z | |
| Neonatal / Infant Silicone Donut Masks | 17 | Zen-O lite™ Portable Oxygen Concentrator | 12 |
| Neonatal Ear Covers, Preemie Muffs | 14 | Zen-O™ Portable Oxygen Concentrator | 11 |







TERMS AND CONDITIONS OF SALE

Parker Healthcare (ABN 43 006 420 364) hereafter known as "the Company".

- 1. All goods purchased in Australia, with the exception of Furniture and orders under \$200.00 are delivered Free Into Store (F.I.S.) (unless otherwise agreed in writing). All goods purchased outside Australia, are supplied at ex-factory prices (unless otherwise agreed in writing).
- 2. All products delivered to On-sellers and Distributors are supplied at ex-Factory Melbourne prices (unless otherwise agreed in writing).
- 3. All risks in the products not delivered F.I.S., shall pass to the "Purchaser" when the product leaves the "Company" store and from that point, regardless of how delivery is affected, any insurance shall be the purchaser's responsibility. Please state Insurance Open Policy number on confirmation order.
- 4. All orders must meet the "Company's" minimum order value of \$200.00. Any order under this limit will attract a "Minimum Order Fee" to the value of the difference between the order and \$200.00.
- 5. Written quotations are firm for 30 days (unless otherwise stated) from date of writing except for material changes in exchange rates or increases in duty on imported products.
- 6. Verbal prices quoted are subject to change without notice and are recommended prices only.
- 7. All orders must be confirmed in writing.
- 8. All claims for short delivery will only be recognized if made in writing 48 hours from date of delivery.
- 9. The "Company" takes no responsibility for incorrect freight methods when freight details are not clearly advised at time of ordering and marked on the confirmation order.
- 10. Payment is strictly net 30 days from date of invoice unless otherwise arranged. It is the option of the "Company" to charge an interest fee of 18.0% p.a. on all amounts outstanding after 30 days.
- 11. Title to the products shall not pass to the "Purchaser" until the full invoice price has been received by the "Company". The "Purchaser" shall hold the products solely as bailee, returnable on demand, pending passing of title.
- 12. A Goods Return Authorisation Number must be requested from the "Company" prior to any item being returned for any reason. This number MUST accompany all documentation and all items returned to the "Company" Stores. The "Company" reserves the right to refuse acceptance of returned items not bearing this authorisation number.
- 13. A 20% minimum restocking charge will apply to all goods returned due to an ordering error by the "Purchaser". Freight on returned goods will be the responsibility of the "Purchaser". Goods held for more than seven (7) days from date of sale will not be accepted for return.
- 14. Goods supplied as a special make or a one off requirement will not be accepted for return under any circumstances.
- 15. An official "company" approval number must be quoted when requesting payment for any repairs carried out by an approved Repairer.
- 16. Disposable products are sold in box quantities.
- 17. Should a standing order be cancelled or reduced in any significant way, the "Company" reserves the right to recover the costs of any 'no charge' items previously delivered as part of the existing standing order arrangement.
- Where a customer changes the requested delivery date of an order, payment for the goods is required to be made in accordance with the approved credit terms of the customer, using the original dispatch date as invoice date.

 Furthermore, if the goods are required to be stored for an extended period, the "Company" reserves the right to charge a holding fee.

Parker Healthcare Ptv Ltd

Head Office & International Department:
16 Redland Drive, Mitcham, VIC., 3132 Australia
Phone: +61 (0)3 9872 0222 Facsimile: +61 (0)3 9873 5911
E-mail: sales@parkerhealth.com.au
Represented in all Australian states



AUSTRALIAN OWNED & OPERATED



Sales and Service Australia / Malaysia / New Zealand

Some products and services may not be available to Malaysian and New Zealand customers due to distribution agreements.

All content subject to change without prior notification.

Parker Healthcare Pty Ltd

16 Redland Drive, Mitcham, VIC., 3132 Australia Phone: +61 3 9872 0222 Facsimile: +61 3 9873 5911 E-mail: sales@parkerhealth.com.au www.parkerhealth.com.au

Parker Healthcare Sdn Bhd

No. 22 Jalan Udang Harimau 2, Medan Niaga, Kepong 51200 Kuala Lumpur Malaysia Phone: +60 3 6263 2848 Facsimilie: +60 3 6262 9848

E-mail: sales@parkerhealth.com.my www.parkerhealth.com.my

Parker Healthcare NZ Ltd

PO Box 101042, NSMC, North Shore City 0745 Auckland NZ Phone: +64 9 415 8524 Facsimile: +64 9 415 9672 E-mail: sales@parkerhealth.co.nz www.parkerhealth.co.nz